


## EZ-Vision® RNA Ladder

*Instant ladder visualization without the use of ethidium bromide*



**EZ-Vision® RNA Ladder.** 7 bands ranging from 1,500 nt - 9,000 nt.

## RNA EZ Vision® Dye as Loading Buffer

AMRESCO's non-mutagenic alternative to ethidium bromide for RNA

### Convenient Solution


Ready-to-use solution for RNA staining in a loading buffer. Recommended for denaturing agarose gels with formaldehyde. Compatible with downstream applications.

### Instant Results

Immediately visualize RNA ladder post-electrophoresis with no staining or destaining required.

### Safer to Use

Premixed solution reduces hazardous chemical handling and completely eliminates the use of ethidium bromide.



**Total RNA was purified from K562 cells with Ribozol™ RNA Extraction Reagent.** Samples were denatured at 65°C for 10 minutes in RNA EZ-Vision® loading buffer containing either 0.003% ethidium bromide or RNA EZ-Vision® dye. The denatured RNA was loaded onto a 2% formaldehyde denaturing agarose gel (Agarose I™, Code: 0710) and resolved for 1.5 hours at 5.1 V/cm. Bands were visualized by UV illumination and image capture was performed with a Syngene GBox-HR Gel Doc System with a SYBR® Green filter.


### Other EZ-Vision Products

- EZ-Vision® In-Gel Solution (for DNA)*
- EZ-Vision® DNA Dye as Loading Buffer*
- EZ-Vision® 1kb DNA Ladder*
- EZ-Vision® 100 bp DNA Ladder*
- EZ-Vision® PCR DNA Ladder*
- Protein EZ-Vision®*

Description	Product Code	Pack Size
<b>RNA EZ-Vision® Dye as Loading Buffer</b>	N717-2X1.25ML	2 X 1.25 ml
<b>EZ-Vision® RNA Ladder</b>	N853-500UL	500 µl

## Formaldehyde-Free RNA Gel Kit

*Instant band visualization  
without the use of  
Ethidium Bromide*



30 25 20 15 10 5 Total RNA (µg)  
← Ethidium Bromide → ← Formaldehyde-Free RNA Dye →

**Comparison of Formaldehyde-Free RNA dye to ethidium bromide staining.** Total RNA was extracted from K562 cells with RiboZol™ RNA Extraction Reagent (N580). Samples were denatured in Formaldehyde-Free RNA Loading Buffer containing either ethidium bromide or Formaldehyde-Free RNA dye and incubated for 10 minutes at 65°C. The samples were loaded on a 2% Formaldehyde-Free RNA Agarose Gel (Agarose I™, Code:0710) and resolved for 1.5 hours at 5.1 V/cm. Image capture was performed with a Syngene GBox-HR Gel Doc System with a SYBR® Green filter.

## RNA Convenience Products from AMRESCO

Description	Product Code	Pack Size
<b>ELECTROPHORESIS</b>		
<b>Formaldehyde-Free RNA Gel Kit</b> Includes gel buffer, running buffer and sample loading buffer. Sample loading buffer includes both an RNA denaturing agent and a non-mutagenic fluorescent dye. Bright green bands are immediately visible upon illumination with UV light without the need for post-run staining or destaining.	N726-KIT	1 Kit
<b>Agarose I™</b> All purpose agarose suitable for analytical and preparative applications; nuclease & protease-free.	0710-500G	500 g
<b>ISOLATION &amp; PURIFICATION</b>		
<b>Phenol-Free Total RNA Purification Kit</b> Purify high quality RNA in just 20 minutes.	N788-KIT	50 Preps
<b>RiboZol™ RNA Extraction Reagent</b> Ready-to-use, single phase extraction reagent.	N580-100ML	100 ml
<b>RiboZol™ Plus RNA Purification Kit</b> Generate high yields of purified, concentrated total RNA with a simple 2 part procedure that combines an initial extraction step with RiboZol™ RNA Extraction Reagent followed by spin column chromatography for additional purification and concentration.	N643-KIT	1 Kit
<b>NucleasEliminator™</b> Removes and deactivates high concentrations of nuclease contamination from glass or plastic surfaces. Available as liquid or single-use wipes.	E891-500ML	500 ml
<b>RiboReserve™ RNA Storage Solution</b> Protects purified RNA from degradation during extended storage. Sodium citrate buffer, pH 6.4 buffer.	N633-12X1ML	12 X 1 ml
<b>TE Buffer, 1X Sterile Solution</b> A biotechnology grade solution for DNA purification and storage.	E112-500ML	500 ml
<b>Water, Sterile, Nuclease-Free</b> Nuclease-free water, specially prepared using DEPC to ensure removal of RNase contamination.	E476-500ML	500 ml