

## ■ Formulations

### DULBECCO'S PHOSPHATE BUFFERED SALINE (D-PBS)

COMPONENTS	Cat.# 114-057 Cat.# 130-057 mg/L	Cat.# 114-059 Cat.# 130-059 mg/L	10X Cat.# 119-068 mg/L
Inorganic Salts			
CaCl <sub>2</sub> •2H <sub>2</sub> O		130.00	
MgCl <sub>2</sub> •6H <sub>2</sub> O		100.00	
KCl	200.00	200.00	2000.00
KH <sub>2</sub> PO <sub>4</sub>	200.00	200.00	2000.00
NaCl	8000.00	8000.00	80000.00
Na <sub>2</sub> HPO <sub>4</sub> •7H <sub>2</sub> O	2160.00	2160.00	21600.00

Reference: Dulbecco, R. and Vogt, M., (1954) *J. Exp. Med.*, **99**:167-182