

Formulations

PHOSPHATE BUFFERED SALINE

COMPONENTS	Cat.# 114-056 mg/L	Cat.# 114-058 Cat.# 130-058 mg/L	10X Cat.# 119-069 Cat.# 130-069 mg/L
Inorganic Salts			
Na ₂ HPO ₄ (anhydrous)	726.00	795.00	7950.00
KH ₂ PO ₄	210.00	144.00	1440.00
NaCl	9000.00	9000.00	90000.00

Reference: Personal Communication - Mr. George Gardiner, NIH Media Unit

DULBECCO'S PHOSPHATE BUFFERED SALINE (D-PBS) WITH TWEEN 20[®]

COMPONENTS	Cat.# 114-154 Cat.# 130-154 mg/L
Inorganic Salts	
KCl	200.00
KH ₂ PO ₄	200.00
NaCl	8000.00
Na ₂ HPO ₄ •7H ₂ O	2160.00
Tween 20 [®]	0.5 ml/L

Tween 20[®] is a registered trademark of ICI Americas, Inc.

Reference: NIH Media Unit