



In the U.S., the chance of a woman being diagnosed with invasive breast cancer in her lifetime is about 1 in 8.  
– 2008 American Cancer Society

## Microflex and The Breast Cancer Research Foundation® *Touching the Lives of Others*


Microflex is proud to be an official Corporate Alliance Partner of The Breast Cancer Research Foundation® by making a donation from every case of *ColorTouch® Pink* sold.

**The mission of the Breast Cancer Research Foundation® is to achieve prevention and a cure for breast cancer in our lifetime** by providing critical funding for innovative clinical and translational research at leading medical centers worldwide, and increasing public awareness about good breast health.

**A minimum of 85% of all funds raised** goes to breast cancer research grants and awareness programs.

BCRF has consecutively received the highest ratings from charity watchdogs for responsible and efficient use of funds.


- Highest rating from Charity Navigator – BCRF outperformed 99% of over 5,000 evaluated charities in efficiently managing its finances
- Named one of the "Top Rated" cancer organizations by The American Institute of Philanthropy – the only breast cancer organization to appear on its list.



**BCRF-funded research has helped save lives and significantly improved the quality of care and survivorship** for tens of thousands of breast cancer patients in the past fifteen years.

BCRF grantees continue to be leaders in uncovering the clues to genetic susceptibility to breast cancer by funding:

- Collection of more than 1300 human genes related to breast cancer, an effort led by Drs. Joan Brugge and Joshua LaBaer at Harvard
- Discovery of new genetic risk markers for breast cancer, including one discovered by Dr. Kenneth Offit at Memorial Sloan-Kettering Cancer Center



BCRF has also been instrumental in the creation and ongoing support of the Translational Breast Cancer Research Consortium (TBCRC), uniting the efforts of the 15 leading breast cancer research centers in the US.

For more information about BCRF, please visit [www.bcrfcure.org](http://www.bcrfcure.org).



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A division of BarrierSafe Solutions International®, Inc.

# Touching the Lives of Others

## ColorTouch® Pink

POWDER-FREE LATEX STANDARD EXAMINATION GLOVES

ColorTouch® Pink latex gloves from Microflex deliver medical-grade barrier protection with consistent comfort and quality in every box. Made from high quality, polymer-coated latex, ColorTouch® Pink latex gloves go on easily and wear comfortably. With relaxing all-day comfort and a wonderful cause attached to its sales, ColorTouch® Pink latex gloves are a worthwhile change from ordinary hand protection.

- Powder-free premium latex
- Polymer-coated for easy donning
- Fully textured for more secure grip
- Relaxed comfort for reduced hand fatigue
- \$1 from each case sold donated to Breast Cancer Research Foundation®

ColorTouch® Pink latex gloves from Microflex are more than great hand protection, they are a product you can be proud to use every day, as your purchases directly support the Breast Cancer Research Foundation® in their search for a cure to breast cancer.



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®



### ColorTouch® Pink Attributes



For current specification information, please visit [www.microflex.com](http://www.microflex.com).

	ASTM D 3578	ColorTouch® Pink
<b>Length (mm/in):</b>	220 / 8.7	245 / 9.6
<b>Cuff Thickness (mils):</b>	n/a	3.9
<b>Palm Thickness (mils):</b>	min. 3.1	5.9
<b>Finger Thickness (mils):</b>	min. 3.1	6.7
<b>Tensile Strength</b>		
Before Aging (MPa):	min. 18	21
After Aging (MPa):	min. 14	18
<b>Elasticity/Elongation:</b>		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600

#### Order #

XS .. CTP-233-XS M ... CTP-233-M  
S .... CTP-233-S L .... CTP-233-L

**Material:** Latex  
**Type:** Non-sterile  
**Color:** Pink

**Cuff:** Beaded  
**Glove Interior:** Powder-free, polymer-coated

**Packaging:** 1000 gloves /case, by weight  
100 gloves /box, by weight, 10 boxes /case

**Usage:** For single use only

**Shape:** Ambidextrous

**Size:** XS - L

**Glove Exterior:** Fully textured

⚠ Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

See for yourself why Microflex® gloves are the first choice for exam glove users nationwide. For more information, or to order a sample, please call: (800) 876-6866. Or visit: [www.microflex.com](http://www.microflex.com)

