

# ColorTouch<sup>®</sup>

LIGHTLY POWDERED LATEX EXAMINATION GLOVES

• STANDARD EXAMINATION GLOVE



## Features & Benefits:

- Premium quality, low modulus latex for greater comfort and reduced hand fatigue
- Lightly powdered with absorbable cornstarch (USP) for easy donning
- Fully textured for a more secure grip
- Meets federal guidelines for maximum powder and maximum water extractable protein content
- Lightly scented with peppermint
- Color-coded sizes

**A Colorful Solution in High Quality Hand Protection.**

**Color Touch<sup>®</sup>** is a refreshing change from ordinary latex gloves with its peppermint scent, and bright, size-identifying colors. Made from premium latex, **Color Touch<sup>®</sup>** goes on easily and wears comfortably. Fully textured for a secure grip in wet and dry conditions, and meets or exceeds current ASTM standards for exam glove quality and performance.

**MICROFLEX<sup>®</sup>**  
THE MOST TRUSTED NAME IN GLOVES<sup>®</sup>

REGISTERED  
**ISO 13485**

# ColorTouch®

## LIGHTLY POWDERED LATEX EXAMINATION GLOVES

### EXAMINATION GLOVE

#### GLOVE PROFILE

	ASTM*	Color Touch®
<b>Length (mm/in):</b>	220 (8.7) / 230 (9.1)	245 / 9.6
<b>Cuff Thickness (mils):</b>	n/a	3.9
<b>Palm Thickness (mils):</b>	min. 3.1	5.9
<b>Finger Thickness (mils):</b>	min. 3.1	6.7
<b>Tensile Strength</b>		
Before Aging (MPa):	min. 18	min. 21
After Aging (MPa):	min. 14	min. 18
<b>Elasticity/Elongation:</b>		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
<b>Reduced Protein Rating:</b>	max. 200 µg/dm <sup>2</sup>	max. 150 µg/g

#### GLOVE ATTRIBUTES

<b>Material:</b>	Premium quality natural rubber latex
<b>Type:</b>	Non-sterile
<b>Shape:</b>	Ambidextrous
<b>Color and Size:</b>	Size coded – Small (pink), Medium (green), Large (blue)
<b>Cuff:</b>	Beaded
<b>Glove Interior:</b>	Lightly powdered with absorbable cornstarch, USP
<b>Glove Exterior:</b>	Textured
<b>Usage:</b>	For single use only
<b>Packaging:</b>	100 gloves per dispenser box, by weight 10 Boxes per case

\* ASTM D 3578



## Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**(800) 876-6866**

Or visit our website at:

**[www.microflex.com](http://www.microflex.com)**

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
<b>CT-133-S</b>	<b>SMALL</b>	Latex, Lightly Powdered, Examination	<b>1000</b>
<b>CT-133-M</b>	<b>MEDIUM</b>	Latex, Lightly Powdered, Examination	<b>1000</b>
<b>CT-133-L</b>	<b>LARGE</b>	Latex, Lightly Powdered, Examination	<b>1000</b>

**CAUTION STATEMENT:** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.