

DIAMOND GRIP PLUS™



POWDER-FREE LATEX EXAM GLOVES

• STANDARD EXAMINATION GLOVE

Features & Benefits:

- Polymer-coated for easy, powder-free donning
- Longer and thicker than typical exam gloves for added protection
- Fully textured for a secure grip in wet or dry conditions
- Premium, low-modulus natural rubber latex for a true fit and reduced hand fatigue
- Reduced protein rating



- Diamond Grip Plus™**
- Powder-Free
 - Textured

Performance and Comfort is a Premium.

Diamond Grip Plus™ is made from soft, premium natural rubber latex for a comfortable fit. Fully textured for enhanced sensitivity and control. Manufactured with integrated polymer technologies for powder-free donning and product consistency. **Diamond Grip Plus™** meets or exceeds current ASTM standards for exam glove quality and performance.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED
5013485

DIAMOND GRIP PLUS™

POWDER-FREE LATEX EXAM GLOVES



STANDARD EXAMINATION GLOVE

GLOVE PROFILE

	ASTM*	Diamond Grip Plus™
Length (mm/in):	220/8.7	248/9.8
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.1
Finger Thickness (mils):	min. 3.1	6.3
Tensile Strength		
Before Aging (MPa):	min. 18	min. 21
After Aging (MPa):	min. 14	min. 18
Elasticity/Elongation:		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
Reduced Protein Rating:	max. 200 µg/dm ²	≤ 50 µg/g***

GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Natural
Sizes:	X-small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Powder-free
Glove Exterior:	Fully textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

* ASTM D 3578

** 2004 Glove Specifications.

*** PROTEIN RATING: This latex glove contains 50 micrograms or less of total water extractable protein per gram.

Diamond Grip Plus™

- Powder-Free
- Textured

Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

(800) 876-6866

Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
DGP-350-XS	EXTRA-SMALL	Latex, Powder-free, Examination	1000
DGP-350-S	SMALL	Latex, Powder-free, Examination	1000
DGP-350-M	MEDIUM	Latex, Powder-free, Examination	1000
DGP-350-L	LARGE	Latex, Powder-free, Examination	1000
DGP-350-XL	EXTRA-LARGE	Latex, Powder-free, Examination	1000

CAUTION STATEMENT: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.