

EVOLUTION ONE®

POWDER-FREE LATEX EXAM GLOVES

• LONGER AND THICKER

Features & Benefits:

- Interior 5-polymer matrix system provides superior anti-sticking properties resulting in highly effective donning in wet, damp or dry conditions
- Fully textured for a secure grip
- Longer beaded cuff for added protection over the wrist and arm
- Premium, low-modulus latex glove for true fit and reduced hand fatigue



- Evolution One®**
- Powder-free
 - Longer cuff

Advanced Polymer Technology. Excellent Barrier Performance.

Evolution One® latex gloves are the standard for exceptional comfort and premium quality in the laboratory setting. This polymer-coated glove is highly durable and goes on easily without powder. Longer and thicker than standard examination gloves for added barrier protection, and fully textured for better gripping in a variety of conditions. **Evolution One®** meets or exceeds current ASTM standards.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED
ISO 13485

EVOLUTION ONE®

POWDER-FREE LATEX EXAM GLOVES



LONGER AND THICKER

GLOVE PROFILE

	ASTM*	Evolution One®**
Length (mm/in):	min. 220/8.7	255/10
Cuff Thickness (mils):	n/a	4.7
Palm Thickness (mils):	min. 3.1	6.7
Finger Thickness (mils):	min. 3.1	8.7
Tensile Strength		
Before Aging (MPa):	min. 18	21
After Aging (MPa):	min. 14	18
Elasticity/Elongation:		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
Reduced Protein Rating:	max. 200 µg/dm ²	≤ 50 µg/g***

GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Natural
Sizes:	X-Small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Polymer-coated, powder-free
Glove Exterior:	Textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

* ASTM D 3578

** 2004 Glove Specifications.

*** PROTEIN RATING: This latex glove contains 50 micrograms or less of total water extractable protein per gram.

- Evolution One®**
- Powder-free
 - Longer cuff

Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

(800) 876-6866

Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
EV-2050-XS	EXTRA-SMALL	Latex, Non-sterile, Examination	1000
EV-2050-S	SMALL	Latex, Non-sterile, Examination	1000
EV-2050-M	MEDIUM	Latex, Non-sterile, Examination	1000
EV-2050-L	LARGE	Latex, Non-sterile, Examination	1000
EV-2050-XL	EXTRA-LARGE	Latex, Non-sterile, Examination	1000

CAUTION STATEMENT: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
2301 ROBB DRIVE • RENO, NV 89523