

FreeForm® EC

POWDER-FREE BLUE NITRILE EXTENDED CUFF EXAMINATION GLOVES

Recommended for use in the following industries:

Laboratory, EMS, Corrections / Law Enforcement,
Industrial Safety, Automotive

With an extended cuff for added protection over the wrist and arm, **FreeForm® EC's** patented* formulation offers comfortable strength in a relaxed fit nitrile examination glove.


Unlike conventional nitrile gloves, which tend to feel stiff and unforgiving, **FreeForm® EC** works with your hands as they move. **FreeForm® EC** is a model example of strength and comfort in a nitrile glove.

FreeForm® EC gloves meet or exceed the emergency medical glove requirements of NFPA 1999, *Standard on Protective Clothing for Emergency Medical Operations*, 2003 edition.**

** Certification on file.

Recommended for use when:

- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and comfort
- Working in environments that demand consistent grip and reliability



NOT MADE FROM NATURAL RUBBER LATEX

FreeForm EC's Attributes



| | ASTM D 6319 | FreeForm® EC |
|---------------------------------|-------------|--------------|
| Length (mm/in): | 220 / 8.7 | 295 / 11.4 |
| Cuff Thickness (mils): | n/a | 4.7 |
| Palm Thickness (mils): | min. 2.0 | 6.3 |
| Finger Thickness (mils): | min. 2.0 | 7.5 |
| Tensile Strength | | |
| Before Aging (MPa): | min. 14 | 18 |
| After Aging (MPa): | min. 14 | 16 |
| Elasticity/Elongation: | | |
| Before Aging (%): | min. 500 | 500 |
| After Aging (%): | min. 400 | 400 |

Order #

S . . . FFE-775-S XL . . FFE-775-XL
M . . . FFE-775-M XXL . . FFE-775-XXL
L FFE-775-L 3XL . . FFE-775-3XL

Material: Nitrile
Type: Non-sterile
Color: Blue
Cuff: Extended, beaded
Glove Interior: Polymer-coated, powder-free
Packaging: 500 gloves per case, by weight, 10 boxes per case
50 gloves per dispenser box, by weight

Usage: For single use only
Shape: Ambidextrous
Size: S - 3XL
Glove Exterior: Textured fingers

⚠ Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

*Patent No. 6,451,893. Patent No. 7,176,260.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. Contact a Microflex® distributor today to purchase quality Microflex® products. For the name of a distributor in your area, call: **(800) 876-6866**. Fax: (800) 876-6632.

