

FreeForm® SE

POWDER-FREE NITRILE EXAM GLOVES

● STANDARD EXAMINATION GLOVE

Unbelievably Nitrile!



FreeForm® SE
Powder-free

FreeForm® SE's patented* formulation offers comfortable strength in a relaxed fit nitrile examination glove.

Unlike conventional nitrile gloves, which tend to feel stiff and unforgiving, FreeForm® SE works with your hands as they move. FreeForm® SE is a model example of strength and comfort.

Features & Benefits:

- Exclusive, patented* formulation
- Not made from Natural Rubber Latex to eliminate the potential for latex protein-mediated allergic sensitization
- Polymer coated to aid in donning ease

*Patent No. 6,451,893

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED
ISO 13485

FreeForm® SE

POWDER-FREE NITRILE EXAM GLOVES



FreeForm® SE
Powder-free

● STANDARD EXAMINATION GLOVE

GLOVE PROFILE

	ASTM*	FreeForm® SE
Length (mm/in):	220/8.7	245/9.6
Cuff Thickness (mils):	n/a	3.5
Palm Thickness (mils):	min. 2.0	3.9
Finger Thickness (mils):	min. 2.0	5.5
Tensile Strength		
Before Aging (MPa):	min. 14	min. 18
After Aging (MPa):	min. 14	min. 16
Elasticity/Elongation:		
Before Aging (%):	min. 500	min. 500
After Aging (%):	min. 400	min. 400

GLOVE ATTRIBUTES

Material:	Nitrile; not made from natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Blue
Sizes:	X-Small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Polymer-coated, powder-free
Glove Exterior:	Textured fingers
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 boxes per case

* ASTM D 6319

Nitrile Examination Gloves

The future of nitrile is now!

See for yourself why **Microflex®** gloves are the first choice of exam glove users. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or call **(800) 876-6866**

for the name of a distributor in your area. Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
FFS-700-XS	EXTRA-SMALL	Nitrile, Non-sterile, Examination	1000
FFS-700-S	SMALL	Nitrile, Non-sterile, Examination	1000
FFS-700-M	MEDIUM	Nitrile, Non-sterile, Examination	1000
FFS-700-L	LARGE	Nitrile, Non-sterile, Examination	1000
FFS-700-XL	EXTRA-LARGE	Nitrile, Non-sterile, Examination	1000

Caution Statement: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
2301 ROBB DRIVE • RENO, NV 89523