



## Safety Data Sheet

Revision Date: 13-March-2025

24 Hour Emergency #: 1-800-255-3924

### 1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

PRODUCT NAME: Magnesium Chloride 2M  
CATALOG NUMBER: 340-034-721  
CHEMICAL NAME: An aqueous solution of Magnesium chloride.

#### Company/Undertaking Identification

Quality Biological Inc.  
7581 Lindbergh Drive  
Gaithersburg, MD 20879  
301-840-9331

### 2. HAZARDS IDENTIFICATION

#### GHS – Classification

The product contains no substances which at their given concentration, are considered to be hazardous to health

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	EC No.	REACH Reg. No.	CAS-No	Weight %	Classification
Magnesium chloride	EEC No. 232-094-6	Not available	7786-30-3	40.7	-

### 4. FIRST AID MEASURES

<b>Skin contact</b>	Immediately wash skin with copious amounts of water followed by washing with soap and copious amounts of water. Remove all contaminated clothing.
<b>Eye contact</b>	Flush eyes with copious amounts of water for at least 15 minutes with eyelids separated and call a physician.
<b>Ingestion</b>	Never give anything by mouth to an unconscious person
<b>Inhalation</b>	If person is unconscious seek emergency medical attention, if person is conscious remove to fresh air and call a physician
<b>Notes to physician</b>	Treat symptomatically.

### 5. FIRE-FIGHTING MEASURES



## 2M Magnesium Chloride

**General:** Wear Self-Contained breathing apparatus in pressure demand, MSHA/NIOSH approved. During a fire, irritating and toxic gases may be generated by thermal decomposition.

**Extinguishing Media:** Use extinguishing media appropriate to surrounding fire conditions..

**Autoignition Temperature:** N/A

**Explosion limits:** Not sensitive.

**Flash Point:** Not flammable.

### 6. ACCIDENTAL RELEASE MEASURES

**Use Personal Protective Equipment:** Including Chemical Splash Goggles, Chemical Resistant Gloves, and appropriate clothing to prevent skin exposure. In addition, a Respiratory protection program that complies with OSHA 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant respirator use.

**Spills and Leaks:** Absorb spill with inert material (e.g. Vermiculite, sand, earth) then place in a suitable container. Clean up spills immediately, observing precautions in the Sections 6 and 8. Provide ventilation.

### 7. HANDLING AND STORAGE

#### Handling Storage

Handle in accordance with good industrial hygiene and safety practice.  
Keep containers tightly closed in a dry, and well-ventilated place with ambient temperature.

### 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

#### Occupational exposure controls

##### Exposure limits

Engineering measures      Ensure adequate ventilation, especially in confined areas

##### Personal protective equipment

Respiratory protection	In case of insufficient ventilation wear suitable respiratory equipment
Hand protection	Protective gloves
Eye protection	Safety glasses with side-shields
Skin and body protection	Lightweight protective clothing.
Hygiene measures	Handle in accordance with good industrial hygiene and safety practice
Environmental exposure Controls	Prevent product from entering drains.

### 9. PHYSICAL AND CHEMICAL PROPERTIES

#### General Information

Form      Liquid



### Important Health Safety and Environmental Information

<b>Boiling point/range</b>	°C No data available	°F No data available
<b>Melting point/range</b>	°C No data available	°F No data available
<b>Flash point</b>	°C No data available	°F No data available
<b>Autoignition temperature</b>	°C No data available	°F No data available
<b>Oxidizing properties</b>	No information available	
<b>Water solubility</b>	No data available	

## 10. STABILITY AND REACTIVITY

<b>Stability</b>	Stable under recommended storage conditions.
<b>Materials to avoid</b>	No information available
<b>Hazardous decomposition products</b>	Hydrogen chloride gas
<b>Polymerization</b>	Hazardous polymerisation does not occur.
<b>Incompatible Products</b>	Strong oxidizing agents.

## 11. TOXICOLOGICAL INFORMATION

### Acute toxicity

Irritating to eye, skin and respiratory system.

### Principle Routes of Exposure/ Potential Health effects

<b>Eyes</b>	No information available
<b>Skin</b>	No information available
<b>Inhalation</b>	No information available
<b>Ingestion</b>	No information available

### Specific effects

<b>Carcinogenic effects</b>	No information available
<b>Mutagenic effects</b>	No information available
<b>Reproductive toxicity</b>	No information available
<b>Sensitization</b>	No information available

### Target Organ Effects

Central nervous system (CNS), Kidneys.

## 12. ECOLOGICAL INFORMATION

### Ecotoxicity

Harmful to aquatic organisms.

Chemical Name	Toxicity to Algae	Toxicity to Fish	Microtox	Daphnia Magna (Water Flea)
Magnesium chloride	EC50 = 2200 mg/L 72 h LC 1970-3880 mg/L	LC50 1970-3880 mg/L Pimephales promelas 96 h LC50= 4210 mg/L Gambusia		EC50 = 1400 mg/L 24 h



		affinis 96 h		
--	--	--------------	--	--

## 13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

## 14. TRANSPORT INFORMATION

### IATA

<b>Proper shipping name</b>	Not classified as dangerous in the meaning of transport regulations
<b>Hazard Class</b>	No information available
<b>Subsidiary Class</b>	No information available
<b>Packing group</b>	No information available
<b>UN-No</b>	No information available

## 15. REGULATORY INFORMATION

### International Inventories

### U.S. Federal Regulations

#### **SARA 313**

This product is not regulated by SARA.

#### **Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)**

This product does not contain HAPs.

### U.S. State Regulations

#### **California Proposition 65**

This product does not contain chemicals listed under Proposition 65

## 16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.



# Quality Biological™

Page 5 of 5

**2M Magnesium Chloride**

**THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT; HOWEVER, IT MAY NOT BE ALL INCLUSIVE. THIS INFORMATION SHOULD BE USED ONLY AS A GUIDE. QUALITY BIOLOGICAL, INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM THE HANDLING OR CONTACT WITH THE ABOVE PRODUCT.**

**End of Safety Data Sheet**