

Formulations

DULBECCO'S PHOSPHATE BUFFERED SALINE (D-PBS)

COMPONENTS	Cat.# 114-057 Cat.# 130-057 mg / L	Cat.# 114-059 Cat.# 130-059 mg / L	10X Cat.# 119-068 mg / L
Inorganic Salts			
CaCl ₂ •2H ₂ O		130.00	
MgCl ₂ •6H ₂ O		100.00	
KCl	200.00	200.00	2000.00
KH ₂ PO ₄	200.00	200.00	2000.00
NaCl	8000.00	8000.00	80000.00
Na ₂ HPO ₄ •7H ₂ O	2160.00	2160.00	21600.00
Reference: Dulbecco, R. and Vogt, M., (1954) J. Exp. Med., 99 :167-182			