

Formulations

LB (LENNOX L) BROTH	
COMPONENTS	Cat.# 340-016 g/ L
NaCl	5.00
Tryptone	10.00
Yeast extract	5.00

LB (LENNOX L) PLATE AGAR 1.5%	
COMPONENTS	Cat.# 340-025 g/ L
NaCl	5.00
Agar	15.00
Tryptone	10.00
Yeast extract	5.00

LB (LENNOX L) TOP AGAR 0.7%	
COMPONENTS	Cat.# 340-004 g/ L
NaCl	5.00
Agar	7.00
Tryptone	10.00
Yeast extract	5.00

M9 CA MEDIUM	
COMPONENTS	Cat.# 340-017 g/ L
CaCl ₂	0.011
KH ₂ PO ₄	3.00
MgSO ₄ •7H ₂ O	0.246
NaCl	0.50
Na ₂ HPO ₄ (anhydrous)	7.00
NH ₄ Cl	1.00
Casamino acids	4.00
D-Glucose	4.00

--