

## Formulations

LB (LENNOX L) BROTH	
COMPONENTS	Cat.# 340-016 g/ L
NaCl	5.00
Tryptone	10.00
Yeast extract	5.00

LB (LENNOX L) PLATE AGAR 1.5%	
COMPONENTS	Cat.# 340-025 g/ L
NaCl	5.00
Agar	15.00
Tryptone	10.00
Yeast extract	5.00

LB (LENNOX L) TOP AGAR 0.7%	
COMPONENTS	Cat.# 340-004 g/ L
NaCl	5.00
Agar	7.00
Tryptone	10.00
Yeast extract	5.00

M9 CA MEDIUM	
COMPONENTS	Cat.# 340-017 g/ L
CaCl <sub>2</sub>	0.011
KH <sub>2</sub> PO <sub>4</sub>	3.00
MgSO <sub>4</sub> •7H <sub>2</sub> O	0.246
NaCl	0.50
Na <sub>2</sub> HPO <sub>4</sub> (anhydrous)	7.00
NH <sub>4</sub> Cl	1.00
Casamino acids	4.00
D-Glucose	4.00

--