

## Formulations

### YT MEDIUM

COMPONENTS	2X Cat.# 340-023 g/ L	Cat.# 340-024 g/ L
NaCl	5.00	5.00
Tryptone	16.00	8.00
Yeast extract	10.00	5.00

**References:**

Miller, J.H., (1972) Experiments in Molecular Genetics, Cold Spring Harbor Laboratory Press

Sambrook, J., Fritsch, E.F., Maniatis, T., (1989) Molecular Cloning, A Laboratory Manual, 2nd. Edition, Cold Spring Harbor Laboratory Press