

Safety Data Sheet

Revision Date: 22-January-2018

24 Hour Emergency #: 1-800-255-3924

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

PRODUCT NAME: Gentamicin 50 mg/ml
CATALOG NUMBER: 120-098-661
CHEMICAL NAME: Aqueous concentrated antibiotic solution of Gentamicin sulfate

Company/Undertaking Identification

Quality Biological Inc.
7581 Lindbergh Drive
Gaithersburg, MD 20879
301-840-9331

2. HAZARDS IDENTIFICATION

GHS – Classification

Classification

Skin Irritant 2 Eye Irritant 2 Respiratory Sensitizer 1 Skin Sensitizer 1 Toxic to Reproduction 2

Signal Word

Danger



Health hazards

Respiratory sensitization Category 1

Skin sensitization Category 1

Hazard Statements

H317 - May cause an allergic skin reaction

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary Statements

Prevention

P261 - Avoid breathing dust/fume/gas/mist/vapors/spray

P272 - Contaminated work clothing should not be allowed out of the workplace

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P284 - In case of inadequate ventilation wear respiratory protection

Response

P302 + P352 - IF ON SKIN: Wash with plenty of soap and water

P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention

P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician

7581 Lindbergh Drive • Gaithersburg, MD 20879

(301) 840-9331 • Toll Free (800) 443-9331 • Fax (301) 840-0743

Email: customerservice@qualitybiological.com • Web Site: www.qualitybiological.com



Quality Biological™

Page 2 of 6

Gentamicin 50 mg/ml

P362 + P364 - Take off all contaminated clothing and wash it before reuse

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	EC No.	CAS-No	Weight %	Classification
Gentamicin, sulfate (salt)	215-778-9	1405-41-0	1-10	-

4. FIRST AID MEASURES

Skin contact	Immediately wash skin with copious amounts of water followed by washing with soap and copious amounts of water. Remove all contaminated clothing.
Eye contact	Flush eyes with copious amounts of water for at least 15 minutes with eyelids separated and call a physician.
Ingestion	Never give anything by mouth to an unconscious person
Inhalation	If person is unconscious seek emergency medical attention, if person is conscious remove to fresh air and call a physician
Notes to physician	Treat symptomatically.

5. FIRE-FIGHTING MEASURES

General: Wear Self-Contained breathing apparatus in pressure demand, MSHA/NIOSH approved. During a fire, irritating and toxic gases may be generated by thermal decomposition.

Extinguishing Media: Dry chemical, CO 2, water spray or alcohol-resistant foam.

Autoignition Temperature: N/A

Explosion limits: Not sensitive.

Flash Point: Not determined.

6. ACCIDENTAL RELEASE MEASURES

Use Personal Protective Equipment: Including Chemical Splash Goggles, Chemical Resistant Gloves, and appropriate clothing to prevent skin exposure. In addition, a Respiratory protection program that complies with OSHA 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant respirator use.

Spills and Leaks: Absorb with sand or vermiculite. Pick up and transfer to properly labeled containers. Ventilate area and wash spill site after material pickup is complete.

7. HANDLING AND STORAGE

Handling	Handle in accordance with good industrial hygiene and safety practice.
Storage	Keep containers tightly closed in a dry, cool and well-ventilated place.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION**Occupational exposure controls**

Exposure limits

Engineering measures Ensure adequate ventilation, especially in confined areas

Personal protective equipment

Respiratory protection In case of insufficient ventilation wear suitable respiratory equipment
Hand protection Protective gloves
Eye protection Safety glasses with side-shields
Skin and body protection Lightweight protective clothing.
Hygiene measures Handle in accordance with good industrial hygiene and safety practice
Environmental exposure
Controls Prevent product from entering drains.

9. PHYSICAL AND CHEMICAL PROPERTIES

General Information

Form Liquid

Important Health Safety and Environmental Information

Boiling point/range	°C No data available	°F No data available
Melting point/range	°C No data available	°F No data available
Flash point	°C No data available	°F No data available
Autoignition temperature	°C No data available	°F No data available
Oxidizing properties	No information available	
Water solubility	No data available	

10. STABILITY AND REACTIVITY

Stability Stable under recommended storage conditions.
Materials to avoid No information available
Hazardous decomposition products Carbon oxides. Nitrogen oxides (NOx). Toxic fumes.
Polymerization Hazardous polymerisation does not occur.
Incompatible Products Strong oxidizing agents.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Causes skin, eye and respiratory tract irritation. Sensitizer.

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Gentamicin, sulfate (salt)	5 g/kg (Rat)		

Chronic Toxicity



Chronic Toxicity

Suspected teratogen.

Principle Routes of Exposure/ Potential Health effects

Eyes	No information available
Skin	No information available
Inhalation	No information available
Ingestion	No information available

Specific effects

Carcinogenic effects	No information available
Mutagenic effects	No information available
Reproductive toxicity	No information available
Sensitization	May cause sensitization of susceptible persons

Target Organ Effects

Kidneys, Ears, Nerves.

12. ECOLOGICAL INFORMATION

Ecotoxicity

Persistence and Degradability

No information available.

Bioaccumulation/Accumulation

No information available.

13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

14. TRANSPORT INFORMATION

IATA

Proper shipping name	Not classified as dangerous in the meaning of transport regulations
Hazard Class	No information available
Subsidiary Class	No information available
Packing group	No information available
UN-No	No information available

15. REGULATORY INFORMATION

International Inventories

TSCA	Does not Comply
EINECS/ELINCS	Complies
ENCS	Complies
IECSC	Complies
KECL	Does not Comply
PICCS	Complies
AICS	Complies

U.S. Federal Regulations**SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

U.S. State Regulations**California Proposition 65**

This product does not contain any Proposition 65 chemicals:

16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT; HOWEVER, IT MAY NOT BE ALL INCLUSIVE. THIS INFORMATION SHOULD BE USED ONLY AS A GUIDE. QUALITY BIOLOGICAL, INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM THE HANDLING OR CONTACT WITH THE ABOVE PRODUCT.

End of Safety Data Sheet