

Safety Data Sheet

Revision Date: 04-August-2016

24 Hour Emergency #: 1-800-255-3924

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

PRODUCT NAME: Protein Running Buffer, 10X (Tris Glycine SDS Running Buffer)
CATALOG NUMBER: 351-132-101
CHEMICAL NAME: An aqueous solution of SDS (Sodium Dodecyl Sulfate, Glycine and Tris base [Tris (Hydroxymethyl) aminomethane]).

Company/Undertaking Identification

Quality Biological Inc.
7581 Lindbergh Drive
Gaithersburg, MD 20879
301-840-9331

2. HAZARDS IDENTIFICATION**GHS – Classification**

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.122)

3. COMPOSITION/INFORMATION ON INGREDIENTS

This material is not considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200). See Chemical Name for composition.

4. FIRST AID MEASURES

Skin contact	Immediately wash skin with copious amounts of water followed by washing with soap and copious amounts of water. Remove all contaminated clothing.
Eye contact	Flush eyes with copious amounts of water for at least 15 minutes with eyelids separated and call a physician.
Ingestion	Never give anything by mouth to an unconscious person
Inhalation	If person is unconscious seek emergency medical attention, if person is conscious remove to fresh air and call a physician
Notes to physician	Treat symptomatically.

5. FIRE-FIGHTING MEASURES

10X Tris Glycin Transfer Buffer

General: Wear Self-Contained breathing apparatus in pressure demand, MSHA/NIOSH approved. During a fire, irritating and toxic gases may be generated by thermal decomposition.

Extinguishing Media: Dry chemical, CO₂, water spray or regular foam.

Autoignition Temperature: N/A

Explosion limits: Not sensitive.

Flash Point: Not determined.

6. ACCIDENTAL RELEASE MEASURES

Use Personal Protective Equipment: Including Chemical Splash Goggles, Chemical Resistant Gloves, and appropriate clothing to prevent skin exposure. In addition, a Respiratory protection program that complies with OSHA 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant respirator use.

Spills and Leaks: Absorb spill with inert material (e.g. Vermiculite, sand, earth) then place in a suitable container. Clean up spills immediately, observing precautions in the Sections 6 and 8. Provide ventilation.

7. HANDLING AND STORAGE

Handling Storage

Handle in accordance with good industrial hygiene and safety practice.
Keep containers tightly closed in a dry, cool and well-ventilated place.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Occupational exposure controls

Exposure limits

Engineering measures Ensure adequate ventilation, especially in confined areas

Personal protective equipment

Respiratory protection	In case of insufficient ventilation wear suitable respiratory equipment
Hand protection	Protective gloves
Eye protection	Safety glasses with side-shields
Skin and body protection	Lightweight protective clothing.
Hygiene measures	Handle in accordance with good industrial hygiene and safety practice
Environmental exposure Controls	Prevent product from entering drains.

9. PHYSICAL AND CHEMICAL PROPERTIES

General Information

Form Liquid

7581 Lindbergh Drive • Gaithersburg, MD 20879

(301) 840-9331 • Toll Free (800) 443-9331 • Fax (301) 840-0743

Email: customerservice@qualitybiological.com • Web Site: www.qualitybiological.com

10X Tris Glycin Transfer Buffer

Important Health Safety and Environmental Information

Boiling point/range	°C No data available	°F No data available
Melting point/range	°C No data available	°F No data available
Flash point	°C No data available	°F No data available
Autoignition temperature	°C No data available	°F No data available
Oxidizing properties	No information available	
Water solubility	No data available	

10. STABILITY AND REACTIVITY

Stability	Stable under recommended storage conditions.
Materials to avoid	No information available
Hazardous decomposition products	Carbon oxides. Ammonia. Nitrogen oxides (NOx).
Polymerization	Hazardous polymerisation does not occur.
Incompatible Products	Strong oxidizing agents.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Irritating to eye, skin and respiratory system.

Principle Routes of Exposure/ Potential Health effects

Eyes	No information available
Skin	No information available
Inhalation	No information available
Ingestion	No information available

Specific effects

Carcinogenic effects	No information available
Mutagenic effects	No information available
Reproductive toxicity	No information available
Sensitization	No information available

Target Organ Effects

None known.

12. ECOLOGICAL INFORMATION

Ecotoxicity

17.92% of the mixture consists of components(s) of unknown hazards to the aquatic environment

13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

14. TRANSPORT INFORMATION

IATA

Not regulated

15. REGULATORY INFORMATION

International Inventories

U.S. Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372..

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

U.S. State Regulations

California Proposition 65

This product does not contain chemicals listed under Proposition 65

16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT; HOWEVER, IT MAY NOT BE ALL INCLUSIVE. THIS INFORMATION SHOULD BE USED ONLY AS A GUIDE. QUALITY BIOLOGICAL, INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM THE HANDLING OR CONTACT WITH THE ABOVE PRODUCT.

End of Safety Data Sheet