

**Safety Data Sheet**

Revision Date: 16-July-2018

24 Hour Emergency #: 1-800-255-3924

**1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**

PRODUCT NAME: LB Broth (Lennox)  
CATALOG NUMBER: 340-004-101, 340-004-131, 340-004-721  
CHEMICAL NAME: Aqueous solution of Sodium chloride, Tryptone, and Yeast extract.

**Company/Undertaking Identification**

Quality Biological Inc.  
7581 Lindbergh Drive  
Gaithersburg, MD 20879  
301-840-9331

**2. HAZARDS IDENTIFICATION****GHS – Classification**

The product contains no substances which at their given concentration, are considered to be hazardous to health

**3. COMPOSITION/INFORMATION ON INGREDIENTS**

Chemical Name	EC No.	REACH Reg. No.	CAS-No	Weight %	Classification
<b>Sodium Chloride</b>	<b>EEC No. 231-598-3</b>	<b>Not available</b>	<b>7647-14-5</b>	<b>1.0</b>	

**4. FIRST AID MEASURES**

<b>Skin contact</b>	Immediately wash skin with copious amounts of water followed by washing with soap and copious amounts of water. Remove all contaminated clothing.
<b>Eye contact</b>	Flush eyes with copious amounts of water for at least 15 minutes with eyelids separated and call a physician.
<b>Ingestion</b>	Never give anything by mouth to an unconscious person
<b>Inhalation</b>	If person is unconscious seek emergency medical attention, if person is conscious remove to fresh air and call a physician
<b>Notes to physician</b>	Treat symptomatically.

**5. FIRE-FIGHTING MEASURES**

## LB Broth (Lennox)

**General:** Wear Self-Contained breathing apparatus in pressure demand, MSHA/NIOSH approved. During a fire, irritating and toxic gases may be generated by thermal decomposition.

**Extinguishing Media:** Dry chemical, CO<sub>2</sub>, water spray or regular foam.

**Autoignition Temperature:** N/A

**Explosion limits:** Not sensitive.

**Flash Point:** Not flammable.

## 6. ACCIDENTAL RELEASE MEASURES

**Use Personal Protective Equipment:** Including Chemical Splash Goggles, Chemical Resistant Gloves, and appropriate clothing to prevent skin exposure. In addition, a Respiratory protection program that complies with OSHA 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant respirator use.

**Spills and Leaks:** Absorb spill with inert material (e.g. Vermiculite, sand, earth) then place in a suitable container. Clean up spills immediately, observing precautions in the Sections 6 and 8. Provide ventilation.

## 7. HANDLING AND STORAGE

### **Handling Storage**

Handle in accordance with good industrial hygiene and safety practice.  
Keep containers tightly closed in a dry, cool and well-ventilated place.

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

### **Occupational exposure controls**

#### **Exposure limits**

**Engineering measures** Ensure adequate ventilation, especially in confined areas

#### **Personal protective equipment**

<b>Respiratory protection</b>	In case of insufficient ventilation wear suitable respiratory equipment
<b>Hand protection</b>	Protective gloves
<b>Eye protection</b>	Safety glasses with side-shields
<b>Skin and body protection</b>	Lightweight protective clothing.
<b>Hygiene measures</b>	Handle in accordance with good industrial hygiene and safety practice
<b>Environmental exposure Controls</b>	Prevent product from entering drains.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### **General Information**

**Form** Liquid

7581 Lindbergh Drive • Gaithersburg, MD 20879

(301) 840-9331 • Toll Free (800) 443-9331 • Fax (301) 840-0743

Email: [customerservice@qualitybiological.com](mailto:customerservice@qualitybiological.com) • Web Site: [www.qualitybiological.com](http://www.qualitybiological.com)

## **Important Health Safety and Environmental Information**

<b>Boiling point/range</b>	°C No data available	°F No data available
<b>Melting point/range</b>	°C No data available	°F No data available
<b>Flash point</b>	°C No data available	°F No data available
<b>Autoignition temperature</b>	°C No data available	°F No data available
<b>Oxidizing properties</b>	No information available	
<b>Water solubility</b>	No data available	

## 10. STABILITY AND REACTIVITY

<b>Stability</b>	Stable under recommended storage conditions.
<b>Materials to avoid</b>	No information available
<b>Hazardous decomposition products</b>	Carbon oxides.
<b>Polymerization</b>	Hazardous polymerisation does not occur.
<b>Incompatible Products</b>	Strong oxidizing agents.

## 11. TOXICOLOGICAL INFORMATION

### **Acute toxicity**

Irritating to eye, skin and respiratory system. May be harmful if swallowed.

### **Principle Routes of Exposure/ Potential Health effects**

<b>Eyes</b>	No information available
<b>Skin</b>	No information available
<b>Inhalation</b>	No information available
<b>Ingestion</b>	No information available

### **Specific effects**

<b>Carcinogenic effects</b>	No information available
<b>Mutagenic effects</b>	No information available
<b>Reproductive toxicity</b>	No information available
<b>Sensitization</b>	No information available

### **Target Organ Effects**

None known.

## 12. ECOLOGICAL INFORMATION

### **Ecotoxicity**

Harmful to aquatic organisms.

Chemical Name	Toxicity to Algae	Toxicity to Fish	Microtox	Daphnia Magna (Water Flea)
Sodium Chloride		LC50 6020-7070 mg/L Pimephales promelas 96 h		EC50 = 1000 mg/L 48 h

		LC50 6420-6700 mg/L Pimephales promelas 96 h LC50 4747-7824 mg/L Oncorhynchus mykiss 96 h LC50= 12946 mg/L Lepomis macrochirus 96 h LC50= 7050 mg/L Pimephales promelas 96 h LC50 5560-6080 mg/L Lepomis macrochirus 96 h		
--	--	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--	--

## 13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

## 14. TRANSPORT INFORMATION

### IATA

<b>Proper shipping name</b>	Not classified as dangerous in the meaning of transport regulations
<b>Hazard Class</b>	No information available
<b>Subsidiary Class</b>	No information available
<b>Packing group</b>	No information available
<b>UN-No</b>	No information available

## 15. REGULATORY INFORMATION

### International Inventories

#### U.S. Federal Regulations

##### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

##### **Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)**

This product does not contain any substances regulated as hazardous air pollutants (HAPS) under Section 112 of the Clean Air Act Amendments of 1990.

#### U.S. State Regulations

##### **California Proposition 65**

This product does not contain chemicals listed under Proposition 65

**WHMIS hazard class:**

D2B Toxic materials



16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.

**THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT; HOWEVER, IT MAY NOT BE ALL INCLUSIVE. THIS INFORMATION SHOULD BE USED ONLY AS A GUIDE. QUALITY BIOLOGICAL, INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM THE HANDLING OR CONTACT WITH THE ABOVE PRODUCT.**

**End of Safety Data Sheet**