

Safety Data Sheet

Revision Date: 04-April-2018

24 Hour Emergency #: 1-800-255-3924

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

PRODUCT NAME: Cell Freezing Medium B (DMSO Freezing Medium 2X)
 CATALOG NUMBER: 118-257-711
 CHEMICAL NAME: Aqueous solution of DMSO, Sodium Chloride, Potassium Chloride, Magnesium Chloride Hexahydrate, Magnesium Sulfate Heptahydrate, Calcium Chloride Anhydrous, Sodium Phosphate Dibasic Heptahydrate, Potassium Phosphate Monobasic, Amino acids, vitamins, Pyruvate, Phenol Red Na, and HEPES.

Company/Undertaking Identification

Quality Biological Inc.
 7581 Lindbergh Drive
 Gaithersburg, MD 20879
 301-840-9331

2. HAZARDS IDENTIFICATION

Hazardous/Non-hazardous Components

The product contains Dimethyl sulfoxide which at their given concentration, are considered to be hazardous to health



Hazard statement(s)

H227 Combustible liquid.

Precautionary statement(s)

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.

P280 Wear protective gloves/ eye protection/ face protection.

P370 + P378 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

P403 + P235 Store in a well-ventilated place. Keep cool.

P501 Dispose of contents/ container to an approved waste disposal plant.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Hazardous components

Component	Classification	Concentration
Dimethyl sulfoxide	Flam. Liq. 4; H227	15%

Principle Routes of Exposure/ Potential Health effects

Eyes	No information available
Skin	No information available
Inhalation	No information available
Ingestion	No information available

Specific effects

Carcinogenic effects	No information available
Mutagenic effects	No information available
Reproductive toxicity	No information available
Sensitization	No information available

Target Organ Effects

No information available

4. FIRST AID MEASURES

Skin contact	Immediately wash skin with copious amounts of water followed by washing with soap and copious amounts of water. Remove all contaminated clothing.
Eye contact	Flush eyes with copious amounts of water for at least 15 minutes with eyelids separated and call a physician.
Ingestion	Never give anything by mouth to an unconscious person
Inhalation	If person is unconscious seek emergency medical attention, if person is conscious remove to fresh air and call a physician
Notes to physician	Treat symptomatically.

5. FIRE-FIGHTING MEASURES

General: Wear Self-Contained breathing apparatus in pressure demand, MSHA/NIOSH approved. During a fire, irritating and toxic gases may be generated by thermal decomposition.

Extinguishing Media: Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Autoignition Temperature: N/A

Explosion limits: N/A

Flash Point: Not Available

6. ACCIDENTAL RELEASE MEASURES

Use Personal Protective Equipment: Including Chemical Splash Goggles, Chemical Resistant Gloves, and appropriate clothing to prevent skin exposure. In addition, a Respiratory protection program that complies with OSHA 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant respirator use.

Spills and Leaks: Absorb spill with inert material (e.g. Vermiculite, sand, earth) then place in a suitable container. Clean up spills immediately, observing precautions in the Sections 6 and 8. Provide ventilation.

7. HANDLING AND STORAGE

Handling Avoid inhalation of vapour or mist. Keep away from sources of ignition - No smoking. Take measures to prevent the build up of electrostatic charge.

Storage Store in a tightly closed container frozen.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Occupational exposure controls

Exposure limits

Engineering measures Ensure adequate ventilation, especially in confined areas

Personal protective equipment

Respiratory protection	In case of insufficient ventilation wear suitable respiratory equipment
Hand protection	Protective gloves
Eye protection	Safety glasses with side-shields
Skin and body protection	Lightweight protective clothing.
Hygiene measures	Handle in accordance with good industrial hygiene and safety practice
Environmental exposure Controls	Prevent product from entering drains.

9. PHYSICAL AND CHEMICAL PROPERTIES

General Information

Form Liquid

Important Health Safety and Environmental Information

Boiling point/range	°C No data available	°F No data available
Melting point/range	°C No data available	°F No data available
Flash point	°C No data available	°F No data available
Autoignition temperature	°C No data available	°F No data available
Oxidizing properties	No information available	
Water solubility	No data available	

10. STABILITY AND REACTIVITY

Reactivity	No data available.
Materials to avoid	No information available
Hazardous decomposition products	Hazardous decomposition products formed under fire conditions. - Carbon oxides, Sulphur oxides
Polymerization	Hazardous polymerisation does not occur.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Principle Routes of Exposure/ Potential Health effects

Eyes	No information available
Skin	No information available
Inhalation	No information available
Ingestion	No information available

Specific effects

Carcinogenic effects	Rat – Oral; Mouse - Oral
Mutagenic effects	No information available
Reproductive toxicity	No information available
Sensitization	No information available

Target Organ Effects

No information available

12. ECOLOGICAL INFORMATION

No Information found.

13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

14. TRANSPORT INFORMATION

IATA

Proper shipping name	Not classified as dangerous in the meaning of transport regulations
Hazard Class	No information available
Subsidiary Class	No information available
Packing group	No information available
UN-No	No information available

15. REGULATORY INFORMATION

International Inventories

U.S. Federal Regulations

SARA 313

This product is not regulated by SARA.

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain HAPs.

U.S. State Regulations

California Proposition 65

This product does not contain chemicals listed under Proposition 65

WHMIS hazard class:

Non-controlled

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR

16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT; HOWEVER, IT MAY NOT BE ALL INCLUSIVE. THIS INFORMATION SHOULD BE USED ONLY AS A GUIDE. QUALITY BIOLOGICAL, INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM THE HANDLING OR CONTACT WITH THE ABOVE PRODUCT.

End of Safety Data Sheet