

Safety Data Sheet

Revision Date: 06-November-2018

24 Hour Emergency #: 1-800-255-3924

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

PRODUCT NAME: 50mg/ml Kanamycin Ready-Made Solution
CATALOG NUMBER: 120-349-731
CHEMICAL NAME: An aqueous solution of Kanamycin Sulfate Salt.

Company/Undertaking Identification

Quality Biological Inc.
7581 Lindbergh Drive
Gaithersburg, MD 20879
301-840-9331

2. HAZARDS IDENTIFICATION**GHS – Classification**

Reproductive toxicity (Category 1B), H360

Signal Word**DANGER!**

Hazard Statements: H360: May damage fertility or the unborn child. Precautionary Statements: P201: Obtain special instructions before use. P202: Do not handle until all safety precautions have been read and understood. P281: Use personal protective equipment as required. P308 + P313: IF exposed or concerned: Get medical advice/ attention. P405: Store locked up. P501: Dispose of contents/ container to an approved waste disposal plant.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Identity: Kanamycin Sulfate
CAS number: 70560-51-9
Molecular Formula: C₁₈H₃₈N₄O₁₅S
Molecular Weight: 582.575 g/mol

4. FIRST AID MEASURES



| | |
|---------------------------|---|
| | 50mg/ml Kanamycin ready-made Solution |
| Skin contact | Immediately wash skin with copious amounts of water followed by washing with soap and copious amounts of water. Remove all contaminated clothing. |
| Eye contact | Flush eyes with copious amounts of water for at least 15 minutes with eyelids separated and call a physician. |
| Ingestion | Never give anything by mouth to an unconscious person |
| Inhalation | If person is unconscious seek emergency medical attention, if person is conscious remove to fresh air and call a physician |
| Notes to physician | Treat symptomatically. |

5. FIRE-FIGHTING MEASURES

General: Wear Self-Contained breathing apparatus in pressure demand, MSHA/NIOSH approved. During a fire, irritating and toxic gases may be generated by thermal decomposition.

Extinguishing Media: Substance is non-combustible, use most readily available; Water Spray, Carbon Dioxide, dry chemical powder or appropriate foam.

Autoignition Temperature: N/A

Explosion limits: N/A

Flash Point: Not Available

6. ACCIDENTAL RELEASE MEASURES

Use Personal Protective Equipment: Including Chemical Splash Goggles, Chemical Resistant Gloves, and appropriate clothing to prevent skin exposure. In addition, a Respiratory protection program that complies with OSHA 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant respirator use.

Spills and Leaks: Absorb spill with inert material (e.g. Vermiculite, sand, earth) then place in a suitable container. Clean up spills immediately, observing precautions in the Sections 6 and 8. Provide ventilation.

7. HANDLING AND STORAGE

Handling Always wear recommended Personal Protective Equipment. No special handling advice required

Storage Store in a tightly closed container at frozen.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Occupational exposure controls

Exposure limits

Engineering measures Ensure adequate ventilation, especially in confined areas

Personal protective equipment

Respiratory protection In case of insufficient ventilation wear suitable respiratory equipment
Hand protection Protective gloves



| | |
|---------------------------------|---|
| Eye protection | Safety glasses with side-shields |
| Skin and body protection | Lightweight protective clothing. |
| Hygiene measures | Handle in accordance with good industrial hygiene and safety practice |
| Environmental exposure | |
| Controls | Prevent product from entering drains. |

9. PHYSICAL AND CHEMICAL PROPERTIES

General Information

| | |
|-------------|--------|
| Form | Liquid |
|-------------|--------|

Important Health Safety and Environmental Information

| | | |
|---------------------------------|--------------------------|----------------------|
| Boiling point/range | °C No data available | °F No data available |
| Melting point/range | °C No data available | °F No data available |
| Flash point | °C No data available | °F No data available |
| Autoignition temperature | °C No data available | °F No data available |
| Oxidizing properties | No information available | |
| Water solubility | soluble | |

10. STABILITY AND REACTIVITY

| | |
|---|--------------------|
| Stability | no data available. |
| Materials to avoid | no data available. |
| Hazardous decomposition products | no data available |
| Polymerization | no data available |

11. TOXICOLOGICAL INFORMATION

Acute toxicity

no data available

Principle Routes of Exposure/ Potential Health effects

| | |
|-------------------|--|
| Eyes | May cause eye irritation with susceptible persons. |
| Skin | May cause skin irritation in susceptible persons. |
| Inhalation | May be harmful by inhalation. |
| Ingestion | May be harmful if swallowed. |

Specific effects

| | |
|------------------------------|---------------|
| Carcinogenic effects | Not available |
| Mutagenic effects | Not available |
| Reproductive toxicity | Not available |
| Sensitization | Not available |

Target Organ Effects

No known effects under normal use conditions.



12. ECOLOGICAL INFORMATION

| | |
|----------------------------|---------------------------|
| Ecotoxicity effects | No information available. |
| Mobility | No information available. |
| Biodegradation | No information available |
| Bioaccumulation | No information available. |

13. DISPOSAL CONSIDERATIONS

Dispose of in accordance with local regulations

14. TRANSPORT INFORMATION

IATA

| | |
|-----------------------------|---|
| Proper shipping name | Not classified as dangerous in the meaning of transport regulations |
| Hazard Class | none |
| Subsidiary Class | none |
| Packing group | none |
| UN-No | none |

15. REGULATORY INFORMATION

International Inventories

U.S. Federal Regulations

SARA 313

This product is not regulated by SARA.

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contains HAPs.

U.S. State Regulations

California Proposition 65

This product does not contain chemicals listed under Proposition 65

16. OTHER INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified

7581 Lindbergh Drive • Gaithersburg, MD 20879

(301) 840-9331 • Toll Free (800) 443-9331 • Fax (301) 840-0743

Email: customerservice@qualitybiological.com • Web Site: www.qualitybiological.com



to handle potentially hazardous chemicals must supervise the use of this material.

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT; HOWEVER, IT MAY NOT BE ALL INCLUSIVE. THIS INFORMATION SHOULD BE USED ONLY AS A GUIDE. QUALITY BIOLOGICAL, INC. SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM THE HANDLING OR CONTACT WITH THE ABOVE PRODUCT.

End of Safety Data Sheet