



CUSTOMER CHANGE NOTIFICATION

Reference: CR2022-43

March 16, 2022

Dear Valued Customer,

This is a formal notification of changes that affect products that may be purchased by your organization.

Change Description

Part Number(s) Affected:	<ul style="list-style-type: none">• 112-013-101 DMEM w/o L-Glutamine. 500ml• 112-014-101 DMEM w/L-Glutamine 500ml
Description of Change:	Tighten the pH range from 6.8-7.4 to 7.0 to 7.4
Reason for Change:	Review of historical data show this range is representative of the process

Implementation and Next Steps

Product manufactured after March 21 2022 will have the new pH range.

We trust this documentation is acceptable for your files. Should you have any questions or concerns regarding this notification, please contact me at curleyj@qualitybiological.com.

Sincerely,

Vice President Quality Assurance