

PRODUCT: LB BROTH (LURIA-BERTANI BROTH, Millers)
CATALOG #: 340-004-131 **Size:** 1000ml
LOT #: 720470

CERTIFICATE OF ANALYSIS

	<u>Specifications</u>	<u>Results</u>
PH:	7.2 ± 0.2	7.250
OSMOLALITY:	410 ± 5% mOsm/Kg of water	415
APPEARANCE:	Light amber, clear to slightly opalescent liquid	pass
STERILITY:	Sterilized using saturated steam under pressure	pass
PERFORMANCE:	Passed E. coli growth performance test.	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	02/2014	
EXPIRATION DATE:	02/2016	



Ronald L. Brown, Ph.D.
Director of Quality Assurance

3/18/14

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.