



MAKING SOLUTIONS WORK

PRODUCT: DEPC TREATED WATER
CATALOG#: 351-068-131 Size: 1000ml
LOT#: 720565

CERTIFICATE OF ANALYSIS

This product is ultra-pure, de-ionized by reverse osmosis and sterile filtered.

	<u>Specifications</u>	<u>Results</u>
RESISTIVITY:	> 15.11 Meg Ohms/cm	pass
APPEARANCE:	Clear, colorless liquid	pass
DNASE:	Not detectable	pass
RNASE:	Not detectable	pass
PROTEASE:	Not detectable	pass
ENDOTOXIN:	<0.06 Eu/ml	pass
RECOMMENDED STORAGE CONDITIONS:	15 – 30°C	
MANUFACTURED DATE:	04/2014	
EXPIRATION DATE:	04/2017	

Ronald L. Brown, Ph.D.
Director of Quality Assurance

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.

7581 Lindbergh Drive • Gaithersburg, MD 20879
(301) 840-9331 • Toll Free (800) 443-9331 • Fax (301) 840-0743
Email: customerservice@qualitybiological.com • Web Site: www.qualitybiological.com