

PRODUCT: M9 CA MEDIUM, pH 7.1
CATALOG #: 340-017-101 **Size:** 500ml
LOT #: 724095

CERTIFICATE OF ANALYSIS

	<u>Range</u>	<u>Result</u>
PH:	7.0-7.2	7.092
OSMOLALITY:	For information only	280 mOsm/kg of water
APPEARANCE:	Clear liquid	pass
STERILITY:	No evidence of bacterial or fungal contamination after 14 days testing on fluid thioglycollate and Tryptic Soy Broth.	pass
RECOMMENDED STORAGE CONDITIONS:	2 – 8°C	
MANUFACTURED DATE:	1/2021	
EXPIRATION DATE:	1/2023	

Michael J. Fisher

02/04/2021

Quality Assurance

Date

FOR RESEARCH USE ONLY
NOT FOR IN-VITRO DIAGNOSTIC USE.